

MDM2 Rabbit anti-Human Polyclonal (pSer395) (PE) Antibody - LS-C225861 - LSBio

CatalogID:	LS-C225861
Target:	MDM2 oncogene, E3 ubiquitin protein ligase
Synonyms:	MDM2 Antibody, ACTFS Antibody, Double minute 2 protein Antibody, Hdm2 Antibody, MDM2 variant FB30 Antibody, HDMX Antibody, Oncoprotein Mdm2 Antibody, p53-binding protein Mdm2 Antibody, MDM2 variant FB28 Antibody
Host	MDM2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MDM2 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	MDM2 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser395 of human MDM2.
Specificity:	Human MDM2
Epitope:	pSer395
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences