

MAP2K1 / MEK1 Rabbit anti-Human Polyclonal (pSer217/Ser221) (FITC) Antibody - LS-C225866 - LSBio	
CatalogID:	LS-C225866
Target:	mitogen-activated protein kinase kinase 1 (MAP2K1)
Synonyms:	MAP2K1 Antibody, ERK activator kinase 1 Antibody, MAPK/ERK kinase 1 Antibody, MAPKK 1 Antibody, MKK1 Antibody, PRKMK1 Antibody, MEK 1 Antibody, MAP kinase kinase 1 Antibody, MAPKK1 Antibody, MEK1 Antibody
Family / Subfamily:	Protein Kinase / MAP2K
Host	MAP2K1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAP2K1 / MEK1 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	MAP2K1 / MEK1 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser217/221 of human MEK1.
Specificity:	Human MAP2K1 / MEK1
Epitope:	pSer217/Ser221
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	