

MAP2K2 / MKK2 / MEK2 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C225884 - LSBio

CatalogID:	LS-C225884
Target:	mitogen-activated protein kinase kinase 2 (MAP2K2)
Synonyms:	MAP2K2 Antibody, ERK activator kinase 2 Antibody, MAP kinase kinase 2 Antibody, MAPK/ERK kinase 2 Antibody, MAPKK 2 Antibody, MKK2 Antibody, MEK 2 Antibody, MEK2 Antibody, MAPKK2 Antibody, PRKMK2 Antibody
Family / Subfamily:	Protein Kinase / MAP2K
Host	MAP2K2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	MAP2K2 / MKK2 / MEK2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MAP2K2 / MKK2 / MEK2 antibody was raised against KLH-conjugated synthetic peptide from internal region of human MAP2K2.
Specificity:	Human MAP2K2 / MKK2 / MEK2
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences