

MINA / MINA53 Rabbit anti-Human Polyclonal (C-Terminus) (AP) Antibody - LS-C225908 - LSBio	
<b>CatalogID:</b>	LS-C225908
<b>Target:</b>	MYC induced nuclear antigen (MINA)
<b>Synonyms:</b>	MINA Antibody, MYC-induced nuclear antigen Antibody, NO52 Antibody, MINA53 Antibody, MYC induced nuclear antigen Antibody, MDIG Antibody, Nucleolar protein 52 Antibody
<b>Host</b>	MINA antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Alkaline Phosphatase (AP)
<b>Immunogen Species:</b>	MINA / MINA53 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	MINA / MINA53 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human MINA.
<b>Specificity:</b>	Human MINA / MINA53
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	