

MINA / MINA53 Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C225915 - LSBio	
CatalogID:	LS-C225915
Target:	MYC induced nuclear antigen (MINA)
Synonyms:	MINA Antibody, MYC-induced nuclear antigen Antibody, NO52 Antibody, MINA53 Antibody, MYC induced nuclear antigen Antibody, MDIG Antibody, Nucleolar protein 52 Antibody
Host	MINA antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	MINA / MINA53 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MINA / MINA53 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human MINA.
Specificity:	Human MINA / MINA53
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	