

EP300 / p300 Rabbit anti-Human Polyclonal (pSer89) (FITC) Antibody - LS-C226060 - LSBio

CatalogID:	LS-C226060
Target:	E1A binding protein p300 (EP300)
Synonyms:	EP300 Antibody, E1A binding protein p300 Antibody, Histone acetyltransferase p300 Antibody, KAT3B Antibody, p300 HAT Antibody, E1A-associated protein p300 Antibody, E1A-binding protein, 300kD Antibody, p300 Antibody, RSTS2 Antibody
Host	EP300 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	EP300 / p300 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	EP300 / p300 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser89 of human EP300.
Specificity:	Human EP300 / p300
Epitope:	pSer89
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences