

EP300 / p300 Rabbit anti-Human Polyclonal (pSer1834) (FITC) Antibody - LS-C226067 - LSBio	
<b>CatalogID:</b>	LS-C226067
<b>Target:</b>	E1A binding protein p300 (EP300)
<b>Synonyms:</b>	EP300 Antibody, E1A binding protein p300 Antibody, Histone acetyltransferase p300 Antibody, KAT3B Antibody, p300 HAT Antibody, E1A-associated protein p300 Antibody, E1A-binding protein, 300kD Antibody, p300 Antibody, RSTS2 Antibody
<b>Host</b>	EP300 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	EP300 / p300 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	EP300 / p300 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1834 of human P300.
<b>Specificity:</b>	Human EP300 / p300
<b>Epitope:</b>	pSer1834
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Uses:</b>	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b> Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	