

EP300 / p300 Rabbit anti-Human Polyclonal (pSer1834) (PE) Antibody - LS-C226073 - LSBio	
CatalogID:	LS-C226073
Target:	E1A binding protein p300 (EP300)
Synonyms:	EP300 Antibody, E1A binding protein p300 Antibody, Histone acetyltransferase p300 Antibody, KAT3B Antibody, p300 HAT Antibody, E1A-associated protein p300 Antibody, E1A-binding protein, 300kD Antibody, p300 Antibody, RSTS2 Antibody
Host	EP300 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	EP300 / p300 antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	EP300 / p300 antibody was raised against kLH-conjugated synthetic phosphopeptide around Ser1834 of human P300.
Specificity:	Human EP300 / p300
Epitope:	pSer1834
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	