

PIK3CD / PI3K Rabbit anti-Human Polyclonal (pTyr524) (PE) Antibody - LS-C226130 - LSBio	
CatalogID:	LS-C226130
Target:	phosphatidylinositol-4,5-bisphosphate 3-kinase, catalytic subunit delta (PIK3CD)
Synonyms:	PIK3CD Antibody, p110D Antibody, PI3-kinase subunit delta Antibody, PI3K Antibody, Phosphoinositide-3-kinase C Antibody, PI3K-delta Antibody, PI3Kdelta Antibody, p110DELTA Antibody, PI3-kinase p110 subunit delta Antibody, PtdIns-3-kinase subunit delta Antibody
Family / Subfamily:	Non-Protein Kinase / not assigned-Non-Protein Kinase
Host	PIK3CD antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	PIK3CD / PI3K antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	PIK3CD / PI3K antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr524 of human PIK3CD.
Specificity:	Human PIK3CD / PI3K.
Epitope:	pTyr524
Reactivity:	Human
Purification:	Immunoaffinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	