

Sprouty 4 / SPRY4 Rabbit anti-Human Polyclonal (pTyr75) (FITC) Antibody - LS-C226400 - LSBio	
<b>CatalogID:</b>	LS-C226400
<b>Target:</b>	sprouty homolog 4 (Drosophila) (SPRY4)
<b>Synonyms:</b>	SPRY4 Antibody, Protein sprouty homolog 4 Antibody, Sprouty 4 Antibody, Sprouty homolog 4 (Drosophila) Antibody, Sprouty 4 Antibody, Spry-4 Antibody
<b>Host</b>	SPRY4 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	Sprouty 4 / SPRY4 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	Sprouty 4 / SPRY4 antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr75 of human SPRY4.
<b>Specificity:</b>	Human Sprouty 4 / SPRY4.
<b>Epitope:</b>	pTyr75
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	