

TNIK Rabbit anti-Human Polyclonal (N-Terminus) (PE) Antibody - LS-C226528 - LSBio

CatalogID:	LS-C226528
Target:	TRAF2 and NCK interacting kinase (TNIK)
Synonyms:	TNIK Antibody, ZC2 Antibody, KIAA0551 Antibody
Family / Subfamily:	Protein Kinase / MSN
Host	TNIK antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TNIK antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TNIK antibody was raised against kLH-conjugated synthetic peptide from the N-terminal region of human TNIK.
Specificity:	Human TNIK.
Epitope:	N-Terminus
Reactivity:	Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences