

TNIK Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C226533 - LSBio	
<b>CatalogID:</b>	LS-C226533
<b>Target:</b>	TRAF2 and NCK interacting kinase (TNIK)
<b>Synonyms:</b>	TNIK Antibody, ZC2 Antibody, KIAA0551 Antibody
<b>Family / Subfamily:</b>	Protein Kinase / MSN
<b>Host</b>	TNIK antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	TNIK antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	TNIK antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human TNIK.
<b>Specificity:</b>	Human TNIK.
<b>Epitope:</b>	C-Terminus
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Uses:</b>	IHC, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	