

USP25 Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C226664 - LSBio	
CatalogID:	LS-C226664
Target:	ubiquitin specific peptidase 25 (USP25)
Synonyms:	USP25 Antibody, Deubiquitinating enzyme 25 Antibody, Ubiquitin thiolesterase 25 Antibody, Ubiquitin thioesterase 25 Antibody, USP on chromosome 21 Antibody, Ubiquitin specific protease 25 Antibody
Family / Subfamily:	Protease / Cysteine C19
Host	USP25 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	USP25 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	USP25 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human USP25.
Specificity:	Human USP25.
Epitope:	C-Terminus
Reactivity:	Human, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	