

LIN28A / LIN28 Rabbit anti-Human Polyclonal (Ser134) (FITC) Antibody - LS-C226835 - LSBio

CatalogID:	LS-C226835
Target:	lin-28 homolog A (C. elegans) (LIN28A)
Synonyms:	LIN28A Antibody, CSDD1 Antibody, Lin-28 homolog (C. elegans) Antibody, LIN-28 Antibody, Lin-28A Antibody, LIN28 Antibody, Protein lin-28 homolog A Antibody, Lin-28 homolog A (C. elegans) Antibody, ZCCHC1 Antibody, RNA-binding protein LIN-28 Antibody
Host	LIN28A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	LIN28A / LIN28 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	LIN28A / LIN28 antibody was raised against kLH-conjugated synthetic peptide surrounding Ser134 of human LIN28.
Specificity:	Human LIN28A / LIN28.
Epitope:	Ser134
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	ICC, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences