

PCDHAC2 Rabbit anti-Human Polyclonal (aa623-652) (FITC) Antibody - LS-C226953 - LSBio	
CatalogID:	LS-C226953
Target:	protocadherin alpha subfamily C, 2 (PCDHAC2)
Synonyms:	PCDHAC2 Antibody, PCDH-ALPHA-C2 Antibody, Protocadherin alpha-C2 Antibody
Family / Subfamily:	Cadherin / not assigned-Cadherin
Host	PCDHAC2 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PCDHAC2 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PCDHAC2 antibody was raised against kLH-conjugated synthetic peptide (aa623-652) of human PCDHAC2.
Specificity:	Human PCDHAC2
Epitope:	aa623-652
Reactivity:	Human
Purification:	Protein A purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	