

HSP90B1 / GP96 / GRP94 Rat anti-Chicken Monoclonal (RPE) (9G10) Antibody - LS-C227176 - LSBio	
CatalogID:	LS-C227176
Target:	heat shock protein 90kDa beta (Grp94), member 1 (HSP90B1)
Synonyms:	HSP90B1 Antibody, ECGP Antibody, Endoplasmic Antibody, gp96 homolog Antibody, gp96 Antibody, GRP-94 Antibody, GRP94 Antibody, Tumor rejection antigen 1 Antibody, TRA1 Antibody
Host	HSP90B1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	9G10
Conjugations:	R. Phycoerythrin (RPE)
Immunogen Species:	HSP90B1 / GP96 / GRP94 antibody was raised against Chicken
Antigen Type:	Purified protein
Immunogen:	HSP90B1 / GP96 / GRP94 antibody was raised against purified Grp94 isolated from chicken oviducts.
Specificity:	Detects a 98kD protein corresponding to the molecular mass of Grp94 on SDS PAGE immunoblots. Does not detect human Hsp90, Grp74, or GrpE from E. coli.
Reactivity:	Chicken, Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Horse, Pig, Rabbit, Sheep, Xenopus
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	