

HSP70 + HSC70 Mouse anti-Chicken Monoclonal (PerCP) (BB70) Antibody - LS-C227185 - LSBio	
CatalogID:	LS-C227185
Target:	Heat Shock Protein HSP70 + HSC70
Host	Heat Shock Protein HSP70 + HSC70 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	BB70
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSP70 + HSC70 antibody was raised against Chicken
Antigen Type:	Native protein
Immunogen:	HSP70 + HSC70 antibody was raised against chicken Hsp70/Hsp90 complex.
Specificity:	Detects ~72 and ~73kD proteins corresponding to the molecular mass of inducible Hsp and Hsc on SDS PAGE immunoblots.
Reactivity:	Chicken, Human, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Rabbit, Sheep, Xenopus, Fish, Drosophila, Yeast, Beluga
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	