

HSPD1 / HSP60 Mouse anti-Human Monoclonal (aa383-447) (PerCP) (LK1) Antibody - LS-C227223 - LSBio

CatalogID:	LS-C227223
Target:	heat shock 60kDa protein 1 (chaperonin) (HSPD1)
Synonyms:	HSPD1 Antibody, 60 kDa chaperonin Antibody, CPN60 Antibody, Heat shock protein 65 Antibody, Heat shock protein 60 Antibody, HSP-60 Antibody, HSP60 Antibody, HuCHA60 Antibody, HLD4 Antibody, HSP65 Antibody, p60 lymphocyte protein Antibody, SPG13 Antibody, Chaperonin 60 Antibody, GROEL Antibody
Host	HSPD1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	LK1
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPD1 / HSP60 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	HSPD1 / HSP60 antibody was raised against recombinant human Hsp60, epitope aa383-447.
Specificity:	Detects an ~60kD protein corresponding to the molecular mass of Hsp60 on SDS PAGE immunoblots.
Epitope:	aa383-447
Reactivity:	Human, Monkey, Mouse, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Sheep, Chicken, Xenopus, Drosophila
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences