

CD74 Mouse anti-Human Monoclonal (PerCP) (PIN.1) Antibody - LS-C227255 - LSBio

CatalogID:	LS-C227255
Target:	CD74 molecule, major histocompatibility complex, class II invariant chain
Synonyms:	CD74 Antibody, CD74 antigen Antibody, DHLAG Antibody, HLADG Antibody, HLA-DR-gamma Antibody, Ia-associated invariant chain Antibody, p33 Antibody, MHC HLA-DR gamma chain Antibody, Ia-GAMMA Antibody
Host	CD74 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	PIN.1
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	CD74 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CD74 antibody was raised against human CD74 invariant chain synthetic peptide.
Specificity:	Detects ~33-35kD protein doublet corresponding to the molecular mass of the p33 and p35 forms of human CD74 on SDS PAGE immunoblots.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences