

HMOX1 / HO-1 Mouse anti-Human Monoclonal (aa1-30) (FITC) (1F12-A6) Antibody - LS-C227336 - LSBio	
CatalogID:	LS-C227336
Target:	heme oxygenase (decycling) 1 (HMOX1)
Synonyms:	HMOX1 Antibody, BK286B10 Antibody, HO-1 Antibody, HSP32 Antibody, Heat shock protein, 32-kD Antibody, Heme oxygenase (decycling) 1 Antibody, Heme oxygenase 1 Antibody, HO1 Antibody
Host	HMOX1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	1F12-A6
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HMOX1 / HO-1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HMOX1 / HO-1 antibody was raised against human HO-1 synthetic peptide, amino acids 1-30.
Specificity:	Detects a 32kD protein, corresponding to the molecular mass of HO-1 on SDS Page immunoblots. Does not cross-react with HO-2.
Epitope:	aa1-30
Reactivity:	Human, Mouse, Bovine
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	