

FKBP5 / FKBP51 Mouse anti-Human Monoclonal (PerCP) (Hi51B) Antibody - LS-C227375 - LSBio	
CatalogID:	LS-C227375
Target:	FK506 binding protein 5 (FKBP5)
Synonyms:	FKBP5 Antibody, 51 kd fk506-binding protein Antibody, 51 kDa FK506-binding protein Antibody, 51 kDa FK506-binding protein 5 Antibody, 51 kDa FKBP Antibody, AIG6 Antibody, Androgen-regulated protein 6 Antibody, FKBP-51 Antibody, FKBP-5 Antibody, FKBP51 Antibody, FF1 antigen Antibody, FK506 binding protein 5 Antibody, p54 Antibody, Ptg-10 Antibody, Rotamase Antibody, T-cell FK506-binding protein Antibody, PPlase Antibody, PPlase FKBP5 Antibody, FK506-binding protein 5 Antibody, FKBP54 Antibody, HSP90-binding immunophilin Antibody
Host	FKBP5 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG
Clone Name:	Hi51B
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	FKBP5 / FKBP51 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FKBP5 / FKBP51 antibody was raised against synthetic peptide corresponding to the residues of human FKBP51.
Specificity:	Detects an ~51kD protein representing FKBP51 in cell lysate. Also detects FKBP51 in whole tissue extracts from rat kidney and rat and mouse testes.
Reactivity:	Human, Mouse, Rat, Dog, Hamster, Rabbit
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014 © 2014 LifeSpan BioSciences	