

FKBP4 / FKBP52 Mouse anti-Human Monoclonal (PerCP) (Hi52C) Antibody - LS-C227389 - LSBio	
CatalogID:	LS-C227389
Target:	FK506 binding protein 4, 59kDa (FKBP4)
Synonyms:	FKBP4 Antibody, 59 kDa immunophilin Antibody, 51 kDa FK506-binding protein Antibody, 52 kDa FK506-binding protein Antibody, 52 kDa FKBP Antibody, FKBP-52 Antibody, FKBP51 Antibody, HBI Antibody, Immunophilin FKBP52 Antibody, FK506 binding protein 4, 59kDa Antibody, FK506-binding protein 4 Antibody, HSP56 Antibody, p52 Antibody, p59 Antibody, Rotamase Antibody, p59 protein Antibody, PPIase Antibody, PPIase FKBP4 Antibody, FKBP-4 Antibody, FKBP52 Antibody, FKBP59 Antibody, HSP binding immunophilin Antibody, HSP-binding immunophilin Antibody
Host	FKBP4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	Hi52C
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	FKBP4 / FKBP52 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	FKBP4 / FKBP52 antibody was raised against synthetic peptide corresponding to the residues of human FKBP52.
Specificity:	Detects an ~52kD protein representing FKBP52 in cell lysate. Also detects FKBP52 in whole tissue extracts from rat kidney and rat and mouse testes. Heavy chain migrates close to FKBP52 on SDS PAGE.
Reactivity:	Human, Mouse, Rat, Dog, Hamster
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC - Paraffin, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	