

LAMP2 / CD107b Rat anti-Mouse Monoclonal (PerCP) (GL2A7) Antibody - LS-C227407 - LSBio	
<b>CatalogID:</b>	LS-C227407
<b>Target:</b>	lysosomal-associated membrane protein 2 (LAMP2)
<b>Synonyms:</b>	LAMP2 Antibody, CD107b Antibody, CD107b antigen Antibody, LAMPB Antibody, LGP110 Antibody, LAMP-2 Antibody
<b>Family / Subfamily:</b>	Cancer
<b>Host</b>	LAMP2 antibody was produced in Rat
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	GL2A7
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	LAMP2 / CD107b antibody was raised against Mouse
<b>Antigen Type:</b>	Purified protein
<b>Immunogen:</b>	LAMP2 / CD107b antibody was raised against purified preparation of mouse liver lysosomal membranes.
<b>Specificity:</b>	Recommended for use in ICC. The antibody will label the presumptive lysosomes and late endosomes in cells that have been permeabilized with saponin.
<b>Reactivity:</b>	Mouse, Rabbit
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ICC, Immunofluorescence, Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	