

MAPK14 / p38 Mouse anti-Human Monoclonal (PerCP) (9F12) Antibody - LS-C227426 - LSBio

CatalogID:	LS-C227426
Target:	mitogen-activated protein kinase 14 (MAPK14)
Synonyms:	MAPK14 Antibody, CSAID-binding protein Antibody, CSBP2 Antibody, CSPB1 Antibody, CSBP1 Antibody, Csaids binding protein Antibody, CSBP Antibody, EXIP Antibody, MAP kinase 14 Antibody, MAP kinase p38 alpha Antibody, MAX-interacting protein 2 Antibody, MAPK 14 Antibody, p38 alpha Antibody, PRKM15 Antibody, PRKM14 Antibody, RK Antibody, SAPK2A Antibody, MAP kinase MXI2 Antibody, Mxi2 Antibody, p38ALPHA Antibody, p38 Antibody, p38 MAP kinase Antibody, p38alpha Exip Antibody
Family / Subfamily:	Protein Kinase / MAPK
Host	MAPK14 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	9F12
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	MAPK14 / p38 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	MAPK14 / p38 antibody was raised against full length recombinant protein expressed in E. coli cells.
Specificity:	Detects a ~38kD protein corresponding to the molecular mass of p38 alpha MAPK on SDS PAGE immunoblots.
Reactivity:	Human, Rat
Purification:	Protein G purified
Presentation:	PBS, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences