

TRAP1 / HSP75 Mouse anti-Human Monoclonal (PerCP) (Trap1-6) Antibody - LS-C227522 - LSBio	
CatalogID:	LS-C227522
Target:	TNF receptor-associated protein 1 (TRAP1)
Synonyms:	TRAP1 Antibody, HSP 75 Antibody, HSP75 Antibody, HSP90L Antibody, TNFR-associated protein 1 Antibody, TRAP-1 Antibody
Host	TRAP1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a,k
Clone Name:	Trap1-6
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	TRAP1 / HSP75 antibody was raised against Human
Antigen Type:	Recombinant protein
Immunogen:	TRAP1 / HSP75 antibody was raised against purified recombinant TRAP1.
Specificity:	Detects Trap1. Detects ~75kD.
Reactivity:	Human
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	