

DNMT / DNMT1 Mouse anti-Human Monoclonal (aa637-650) (PerCP) (60B1220.1) Antibody - LS-C227539 - LSBio	
CatalogID:	LS-C227539
Target:	DNA (cytosine-5-)-methyltransferase 1 (DNMT1)
Synonyms:	DNMT1 Antibody, AIM Antibody, CXXC finger protein 9 Antibody, CXXC9 Antibody, DNA methyltransferase 1 Antibody, DNMT Antibody, DNA methyltransferase Hsa1 Antibody, DNA MTase Hsa1 Antibody, HSN1E Antibody, M.Hsa1 Antibody, MCMT Antibody
Family / Subfamily:	C5-methyltransferase
Host	DNMT1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1,k
Clone Name:	60B1220.1
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	DNMT / DNMT1 antibody was raised against Human
Immunogen:	DNMT / DNMT1 antibody was raised against synthetic peptide corresponding to amino acids 637-650 of human DNMT1 (Gene accession number NP_001370).
Specificity:	It will cross-react with mouse DNMT1.
Epitope:	aa637-650
Reactivity:	Human, Mouse, Zebrafish
Purification:	Protein G purified
Presentation:	PBS, 0.05% BSA, 0.05% sodium azide.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation, Chromatin Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	