

HCN1 Mouse anti-Rat Monoclonal (aa778-910) (PerCP) (S70-28) Antibody - LS-C227603 - LSBio	
<b>CatalogID:</b>	LS-C227603
<b>Target:</b>	hyperpolarization activated cyclic nucleotide-gated potassium channel 1 (HCN1)
<b>Synonyms:</b>	HCN1 Antibody, BCNG1 Antibody, BCNG-1 Antibody, HAC2 Antibody, HAC-2 Antibody
<b>Family / Subfamily:</b>	Ion Channel / Potassium channel - cyclic nucleotide gated
<b>Host</b>	HCN1 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG1
<b>Clone Name:</b>	S70-28
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	HCN1 antibody was raised against Rat
<b>Antigen Type:</b>	Fusion protein
<b>Immunogen:</b>	HCN1 antibody was raised against fusion protein amino acids 778-910 (C terminus) of rat HCN1, accession number Q9JKBO.
<b>Specificity:</b>	Detects ~100kD. No cross-reactivity against HCN2.
<b>Epitope:</b>	aa778-910
<b>Reactivity:</b>	Rat, Human, Mouse
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	