

TRPV3 Mouse anti-Rat Monoclonal (aa458-474) (APC) (S15-4) Antibody - LS-C227689 - LSBio	
CatalogID:	LS-C227689
Target:	transient receptor potential cation channel, subfamily V, member 3 (TRPV3)
Synonyms:	TRPV3 Antibody, OLMS Antibody, VRL-3 Antibody, VRL3 Antibody, Vanilloid receptor-like 3 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPV3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	S15-4
Conjugations:	Allophycocyanin (APC)
Immunogen Species:	TRPV3 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	TRPV3 antibody was raised against synthetic peptide amino acids 458-474 (C-terminus) of rat TrpV3, accession number NP_001020928.
Specificity:	Detects ~70kD.
Epitope:	aa458-474
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	