

HCN4 Mouse anti-Rat Monoclonal (aa1019-1198) (PerCP) (S114-10) Antibody - LS-C227701 - LSBio

CatalogID:	LS-C227701
Target:	hyperpolarization activated cyclic nucleotide-gated potassium channel 4 (HCN4)
Synonyms:	HCN4 Antibody, BCNG-3 Antibody, SSS2 Antibody
Family / Subfamily:	Ion Channel / Potassium channel - cyclic nucleotide gated
Host	HCN4 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S114-10
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HCN4 antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	HCN4 antibody was raised against fusion protein amino acids 1019-1198 (c-terminus) of rat HCN4, accession number Q9JKA7.
Specificity:	Detects ~130kD. No cross-reactivity against other HCNs.
Epitope:	aa1019-1198
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences