

TRPV3 Mouse anti-Rat Monoclonal (aa458-474) (AP) (S15-39) Antibody - LS-C227740 - LSBio

CatalogID:	LS-C227740
Target:	transient receptor potential cation channel, subfamily V, member 3 (TRPV3)
Synonyms:	TRPV3 Antibody, OLMS Antibody, VRL-3 Antibody, VRL3 Antibody, Vanilloid receptor-like 3 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPV3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	S15-39
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	TRPV3 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	TRPV3 antibody was raised against synthetic peptide amino acids 458-474 (cytoplasmic C-terminus) of rat TrpV3 accession number NP_001020928).
Specificity:	Detects ~70kD.
Epitope:	aa458-474
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences