

TRPV3 Mouse anti-Rat Monoclonal (aa458-474) (PerCP) (S15-39) Antibody - LS-C227746 - LSBio	
<b>CatalogID:</b>	LS-C227746
<b>Target:</b>	transient receptor potential cation channel, subfamily V, member 3 (TRPV3)
<b>Synonyms:</b>	TRPV3 Antibody, OLMS Antibody, VRL-3 Antibody, VRL3 Antibody, Vanilloid receptor-like 3 Antibody
<b>Family / Subfamily:</b>	Ion Channel / Calcium channel - TRP
<b>Host</b>	TRPV3 antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgG2a
<b>Clone Name:</b>	S15-39
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	TRPV3 antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	TRPV3 antibody was raised against synthetic peptide amino acids 458-474 (cytoplasmic C-terminus) of rat TrpV3 accession number NP_001020928).
<b>Specificity:</b>	Detects ~70kD.
<b>Epitope:</b>	aa458-474
<b>Reactivity:</b>	Rat, Human, Mouse
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ICC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	