

GRIN2B / NMDAR2B / NR2B Mouse anti-Rat Monoclonal (aa20-271) (PerCP) (S59-36) Antibody - LS-C227786 - LSBio	
CatalogID:	LS-C227786
Target:	glutamate receptor, ionotropic, N-methyl D-aspartate 2B (GRIN2B)
Synonyms:	GRIN2B Antibody, GluN2B Antibody, GluR epsilon 2 Antibody, NMDA receptor 2B Antibody, NMDAR2B Antibody, NR2B Antibody, NR3 Antibody, MRD6 Antibody, HNR3 Antibody
Family / Subfamily:	Ion Channel / Glutamate receptor - ionotropic (NMDA receptor)
Host	GRIN2B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	S59-36
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	GRIN2B / NMDAR2B / NR2B antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	GRIN2B / NMDAR2B / NR2B antibody was raised against fusion protein amino acids 20-271 (extracellular N-terminus) of rat NR2B, accession number QOO960.
Specificity:	Detects ~166kD. No cross-reactivity against NR2A.
Epitope:	aa20-271
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	