

TRPC7 Mouse anti-Human Monoclonal (aa845-862) (PerCP) (S64A-36) Antibody - LS-C227809 - LSBio	
CatalogID:	LS-C227809
Target:	transient receptor potential cation channel, subfamily C, member 7 (TRPC7)
Synonyms:	TRPC7 Antibody, KNP3 Antibody, Trrp8 Antibody, Transient receptor protein 7 Antibody, TRP-7 Antibody, TRP7 Antibody, HTRP7 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPC7 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S64A-36
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	TRPC7 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TRPC7 antibody was raised against synthetic peptide amino acids 845-862 of human TRPC7, accession number Q9HCX4.
Specificity:	Detects ~100kD.
Epitope:	aa845-862
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	