

TRPC5 Mouse anti-Human Monoclonal (aa827-845) (PerCP) (S67-15) Antibody - LS-C227827 - LSBio	
CatalogID:	LS-C227827
Target:	transient receptor potential cation channel, subfamily C, member 5 (TRPC5)
Synonyms:	TRPC5 Antibody, CCE2 Antibody, HTRP-5 Antibody, HTRP5 Antibody, Transient receptor protein 5 Antibody, TRP5 Antibody, TRP-5 Antibody, Trrp5 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPC5 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	S67-15
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	TRPC5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TRPC5 antibody was raised against synthetic peptide amino acids 827-845 of human TrpC5 (also known as short transient receptor potential channel 5, and Htrp5.
Specificity:	Detects ~110kD.
Epitope:	aa827-845
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Immunofluorescence, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	