

MNK / ATP7A Mouse anti-Human Monoclonal (aa42-61) (PerCP) (S60-4) Antibody - LS-C228588 - LSBio	
CatalogID:	LS-C228588
Target:	ATPase, Cu ⁺⁺ transporting, alpha polypeptide (ATP7A)
Synonyms:	ATP7A Antibody, Copper pump 1 Antibody, Menkes disease ATPase Antibody, MNK Antibody, MK Antibody, SMAX3 Antibody, Menkes syndrome Antibody, OHS Antibody, Copper-transporting ATPase 1 Antibody, DSMAX Antibody, MC1 Antibody
Family / Subfamily:	Transporter / ATPase - P type, type IB
Host	ATP7A antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	S60-4
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	MNK / ATP7A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	MNK / ATP7A antibody was raised against synthetic peptide amino acids 42-61 (cytoplasmic C-terminus) of human Copper- transporting ATPase 1.
Specificity:	Detects ~180kD in rat brain membrane preparations.
Epitope:	aa42-61
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	