

WC1 / ATP7B Mouse anti-Human Monoclonal (aa3-21) (FITC) (S62-29) Antibody - LS-C228602 - LSBio	
CatalogID:	LS-C228602
Target:	ATPase, Cu++ transporting, beta polypeptide (ATP7B)
Synonyms:	ATP7B Antibody, Copper-transporting ATPase 2 Antibody, WC1 Antibody, Wilson disease Antibody, PWD Antibody, WD Antibody, Copper pump 2 Antibody, WND Antibody
Family / Subfamily:	Transporter / ATPase - P type, type IB
Host	ATP7B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S62-29
Conjugations:	Fluorescein (FITC)
Immunogen Species:	WC1 / ATP7B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	WC1 / ATP7B antibody was raised against synthetic peptide amino acids 3-21 (cytoplasmic N-terminus) of human Copper-transporting ATPase2.
Specificity:	Detects ~160kD in rat brain membrane preparations.
Epitope:	aa3-21
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	