

CIC-3 / CLCN3 Mouse anti-Rat Monoclonal (aa98-115) (PerCP) (S258-5) Antibody - LS-C228703 - LSBio	
CatalogID:	LS-C228703
Target:	chloride channel, voltage-sensitive 3 (CLCN3)
Synonyms:	CLCN3 Antibody, Chloride channel protein 3 Antibody, CIC-3 Antibody, Chloride channel 3 Antibody, Chloride transporter CIC-3 Antibody, CLC3 Antibody
Family / Subfamily:	Ion Channel / Chloride channel - CLC
Host	CLCN3 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S258-5
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	CIC-3 / CLCN3 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	CIC-3 / CLCN3 antibody was raised against synthetic peptide amino acids 98-115 (cytoplasmic N-terminus) of rat Clcn3
Specificity:	Detects ~90kD. Does not cross-react with Clcn4 or Clcn5 (based on KO validation results).
Epitope:	aa98-115
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	