

BCAN / Brevican Mouse anti-Rat Monoclonal (aa219-655) (PerCP) (S294A-6) Antibody - LS-C228920 - LSBio

CatalogID:	LS-C228920
Target:	brevican (BCAN)
Synonyms:	BCAN Antibody, Brevican Antibody, BEHAB Antibody, CSPG7 Antibody, Brevican core protein Antibody
Host	BCAN antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	S294A-6
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	BCAN / Brevican antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	BCAN / Brevican antibody was raised against fusion protein amino acids 219-655 of rat Brevican, accession number P55068.
Specificity:	Detects ~140kD (and smaller due to proteolytic cleavage).
Epitope:	aa219-655
Reactivity:	Rat, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences