

GRIN2A / NMDAR2A / NR2A Mouse anti-Rat Monoclonal (aa75-325) (PerCP) (S327-95) Antibody - LS-C228937 - LSBio	
CatalogID:	LS-C228937
Target:	glutamate receptor, ionotropic, N-methyl D-aspartate 2A (GRIN2A)
Synonyms:	GRIN2A Antibody, EPND Antibody, HNR2A Antibody, GluN2A Antibody, GluR epsilon 1 Antibody, NMDA receptor 2A Antibody, NMDAR2A Antibody, NR2A Antibody, NMDA receptor subtype 2A Antibody
Family / Subfamily:	Ion Channel / Glutamate receptor - ionotropic (NMDA receptor)
Host	GRIN2A antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	S327-95
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	GRIN2A / NMDAR2A / NR2A antibody was raised against Rat
Antigen Type:	Fusion protein
Immunogen:	GRIN2A / NMDAR2A / NR2A antibody was raised against fusion protein amino acids 75-325 (extracellular N-terminus) of rat GluN2A/NR2A, accession number Q00959.
Specificity:	Detects ~170kD. Does not react with NR2B.
Epitope:	aa75-325
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	