

NOTCH1 Mouse anti-Mouse Monoclonal (aa20-216) (PerCP) (S253-32) Antibody - LS-C228955 - LSBio

CatalogID:	LS-C228955
Target:	notch 1 (NOTCH1)
Synonyms:	NOTCH1 Antibody, Notch 1 Antibody, TAN1 Antibody
Host	NOTCH1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Clone Name:	S253-32
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	NOTCH1 antibody was raised against Mouse
Antigen Type:	Fusion protein
Immunogen:	NOTCH1 antibody was raised against fusion protein amino acids 20-216 (extracellular N-terminus, EGF-like domains 1-5) of mouse Notch1, accession number Q01705.
Specificity:	Detects >270kD, ~120kD and small fragments due to proteolysis. Does not cross-react with Notch2 or Notch3.
Epitope:	aa20-216
Reactivity:	Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences