

ATXN1 / Ataxin-1 / SCA1 Mouse anti-Human Monoclonal (aa164-197) (PerCP) (S76-8) Antibody - LS-C229238 - LSBio	
CatalogID:	LS-C229238
Target:	ataxin 1 (ATXN1)
Synonyms:	ATXN1 Antibody, ATX1 Antibody, Ataxin 1 Antibody, D6S504E Antibody, Ataxin-1 Antibody, SCA1 Antibody
Host	ATXN1 antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Clone Name:	S76-8
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	ATXN1 / Ataxin-1 / SCA1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ATXN1 / Ataxin-1 / SCA1 antibody was raised against synthetic peptide amino acids 164-197 (ATTPSQRSQLEAYSTLLANMGSLSQAPGHKVEPP) of mouse Ataxin-1 (Accession No. P54254). Rat: 100% identity (34/34 amino acids identical). Human: 88% identity (30/34 amino acids identical). Percent identity by BLAST analysis: Mouse, Rat (100%).
Specificity:	Detects a 85 kD protein.
Epitope:	aa164-197
Reactivity:	Human, Mouse, Rat
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.1% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	