

TCP1 Rat anti-Mouse Monoclonal (C-Terminus) (PerCP) (91a) Antibody - LS-C229453 - LSBio

CatalogID:	LS-C229453
Target:	t-complex 1 (TCP1)
Synonyms:	TCP1 Antibody, CCTA Antibody, CCT-alpha Antibody, D6S230E Antibody, T-complex 1 Antibody, TCP-1-alpha Antibody, Tailless complex polypeptide 1 Antibody, CCT1 Antibody
Host	TCP1 antibody was produced in Rat
Clonality:	Monoclonal
Isotype:	IgG2a
Clone Name:	91a
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	TCP1 antibody was raised against Mouse
Antigen Type:	Recombinant protein
Immunogen:	TCP1 antibody was raised against recombinant Mouse TCP1 alpha protein fragment (carboxy terminal region).
Specificity:	Detects a 60kD protein corresponding to the molecular weight of TCP1 alpha; also detects a molecular mass band of ~92kD. Cross reactivity with human Hsp60 has been observed with this antibody in immunoblot analysis.
Epitope:	C-Terminus
Reactivity:	Mouse, Human, Monkey, Rat, Bovine, Dog, Guinea pig, Hamster, Pig, Rabbit, Yeast, C. elegans, Fruit fly
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.1% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences