

HSPD1 / HSP60 Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C229470 - LSBio

<b>CatalogID:</b>	LS-C229470
<b>Target:</b>	heat shock 60kDa protein 1 (chaperonin) (HSPD1)
<b>Synonyms:</b>	HSPD1 Antibody, 60 kDa chaperonin Antibody, CPN60 Antibody, Heat shock protein 65 Antibody, Heat shock protein 60 Antibody, HSP-60 Antibody, HSP60 Antibody, HuCHA60 Antibody, HLD4 Antibody, HSP65 Antibody, p60 lymphocyte protein Antibody, SPG13 Antibody, Chaperonin 60 Antibody, GROEL Antibody
<b>Host</b>	HSPD1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	HSPD1 / HSP60 antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	HSPD1 / HSP60 antibody was raised against human HSP60 produced through recombinant DNA methods in E. coli.
<b>Specificity:</b>	Recognizes a ~60kD protein, corresponding to the molecular mass of hsp60 on SDS-PAGE blots.
<b>Reactivity:</b>	Human, Mouse, Rat, Bovine, Dog, Rabbit, Chicken
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, Immunoprecipitation, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences