

KDEL Tetrapeptide Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C229483 - LSBio	
CatalogID:	LS-C229483
Target:	KDEL Tetrapeptide
Host	Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	Human
Antigen Type:	Synthetic peptide
Immunogen:	KDEL containing peptide immunogen
Specificity:	Recognizes KDEL proteins, GRP94, Grp78, PDI and Calreticulin. It may also see ERp57 and ERp29.
Reactivity:	Human, Mouse, Rat
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	200 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	