

HMOX1 / HO-1 Rabbit anti-Human Polyclonal (PerCP) Antibody - LS-C229500 - LSBio	
CatalogID:	LS-C229500
Target:	heme oxygenase (decycling) 1 (HMOX1)
Synonyms:	HMOX1 Antibody, BK286B10 Antibody, HO-1 Antibody, HSP32 Antibody, Heat shock protein, 32-kD Antibody, Heme oxygenase (decycling) 1 Antibody, Heme oxygenase 1 Antibody, HO1 Antibody
Host	HMOX1 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HMOX1 / HO-1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HMOX1 / HO-1 antibody was raised against human heme-oxygenase (HO-1) synthetic multiple antigenic peptide.
Specificity:	Detects an ~32kD protein, corresponding to the molecular mass of HO-1 on SDS Page immunoblots.
Reactivity:	Human, Mouse, Rat, Dog
Purification:	Affinity purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/28/2014 © 2014 LifeSpan BioSciences	