

HSPA5 / GRP78 / BIP Rabbit anti-Rat Polyclonal (PerCP) Antibody - LS-C229617 - LSBio	
CatalogID:	LS-C229617
Target:	heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa) (HSPA5)
Synonyms:	HSPA5 Antibody, BIP Antibody, Heat shock 70 kDa protein 5 Antibody, MIF2 Antibody, GRP-78 Antibody, GRP78 Antibody
Host	HSPA5 antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	HSPA5 / GRP78 / BIP antibody was raised against Rat
Antigen Type:	Recombinant protein
Immunogen:	HSPA5 / GRP78 / BIP antibody was raised against his-tagged rat GRP78.
Specificity:	Detects a 78kD protein corresponding to the molecular mass of Grp78 in SDS PAGE immunoblots.
Reactivity:	Rat, Human, Mouse
Purification:	Protein G purified
Presentation:	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	