

GSK3B / GSK3 Beta Rabbit anti-Human Polyclonal (C-Terminus) (PerCP) Antibody - LS-C229634 - LSBio

CatalogID:	LS-C229634
Target:	glycogen synthase kinase 3 beta (GSK3B)
Synonyms:	GSK3B Antibody, GSK-3 beta Antibody, GSK3 beta Antibody, GSK3beta Antibody, GSK3beta isoform Antibody, TPKI Antibody, GSK 3beta Antibody, Tau protein kinase I Antibody, ZGSK-3beta Antibody
Family / Subfamily:	Protein Kinase / GSK-3
Host	GSK3B antibody was produced in Rabbit
Clonality:	Polyclonal
Conjugations:	Peridinin-chlorophyll-protein (PerCP)
Immunogen Species:	GSK3B / GSK3 Beta antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	GSK3B / GSK3 Beta antibody was raised against synthetic peptide corresponding to the sequence near the C-terminus of human GSK-3Beta.
Specificity:	Detects ~47kD.
Epitope:	C-Terminus
Reactivity:	Human, Mouse, Rat, Bovine
Purification:	Affinity purified
Presentation:	PBS, pH 7.4.
Recommended Storage:	Store at -20°C.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences