

**Ras (H, K, N) Rabbit anti-Human Polyclonal (aa31-43) (PerCP) Antibody - LS-C229651 - LSBio**

<b>CatalogID:</b>	LS-C229651
<b>Target:</b>	Ras (H-, K-, and N-Ras)
<b>Host</b>	Ras (H-, K-, and N-Ras) antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Ras (H, K, N) antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	Ras (H, K, N) antibody was raised against synthetic peptide corresponding to amino acids 31-43 of human RAS. (The immunogen peptide is from the effector binding loop (I2) of the H, K and N-Ras. Sequence is identical yeast, slime mold, fungi, Xenopus, rat, mouse and chicken over residues).
<b>Specificity:</b>	Detects ~21kD.
<b>Epitope:</b>	aa31-43
<b>Reactivity:</b>	Human, Mouse, Rat, Bovine
<b>Purification:</b>	Protein G purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences