

**HSPA5 / GRP78 / BIP Rabbit anti-Human Polyclonal (N-Terminus) (PerCP) Antibody - LS-C229666  
- LSBio**

<b>CatalogID:</b>	LS-C229666
<b>Target:</b>	heat shock 70kDa protein 5 (glucose-regulated protein, 78kDa) (HSPA5)
<b>Synonyms:</b>	HSPA5 Antibody, BIP Antibody, Heat shock 70 kDa protein 5 Antibody, MIF2 Antibody, GRP-78 Antibody, GRP78 Antibody
<b>Host</b>	HSPA5 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	HSPA5 / GRP78 / BIP antibody was raised against Human
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	HSPA5 / GRP78 / BIP antibody was raised against full length human GRP78 (Bip) his tagged at the N terminus.
<b>Specificity:</b>	Detects a protein with a mass of ~78kD corresponding to Grp78.
<b>Epitope:</b>	N-Terminus
<b>Reactivity:</b>	Human, Mouse, Rat, Dog
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences