

HSP70 / Heat Shock Protein 70 Rabbit anti-Plasmodium falciparum Polyclonal (aa365-681) (PerCP) Antibody - LS-C229678 - LSBio	
<b>CatalogID:</b>	LS-C229678
<b>Target:</b>	HSP70 / Heat Shock Protein 70
<b>Host</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>Conjugations:</b>	Peridinin-chlorophyll-protein (PerCP)
<b>Immunogen Species:</b>	Plasmodium falciparum
<b>Antigen Type:</b>	Recombinant protein
<b>Immunogen:</b>	His-tagged and purified PfHsp70, C-terminus (AA 365-681).
<b>Specificity:</b>	This protein A Affinity purified antibody detects an 70kD protein, corresponding to P. falciparum Hsp70 on SDS PAGE immunoblot; in samples from P. falciparum origin.
<b>Epitope:</b>	aa365-681
<b>Reactivity:</b>	Plasmodium falciparum
<b>Non-Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.4, 0.09% sodium azide, 50% glycerol.
<b>Recommended Storage:</b>	Store at -20°C.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Immunofluorescence, Western blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/28/2014	
© 2014 LifeSpan BioSciences	